## **MEGATRON-PNV-125-LF**

## **Nucleus Proton Digital Mixing Valve**

- 1. Nucleus Proton Digital Mixing Valve
  - A. Mixing Valve shall be Digital and comply with National Low Lead Laws @<.25% Lead
  - B. MEGATRON-PNV-125-LF shall:
    - 1. Have 1" inlet/ 1 ¼" outlet connections with integral stop/check valves
    - 2. 2-line, 16-character LED display
    - 3. Integral RTD Sensor
    - 4. Deliver mixed water flow of 95 GPM @ 30 PSI Pressure Drop
    - 5. Maintain temperature with .25 GPM flow from the domestic hot water loop
  - C. MEGATRON-PNV-125-LF shall have all the following operational capabilities:
    - 1. 1X per day shuttle sweep designed to prevent scale buildup on internal mechanical components
    - 2. +/- 2°F water temperature control at times of use and no demand
    - 3. 2°F minimum inlet to outlet water temperature differential
    - 4. Automatic shutoff of hot water upon cold water inlet supply failure
    - 5. Automatic shutoff of cold water upon hot water inlet supply failure
    - 6. Maintain last control position in the event of power failure or be equipped with UPS standby power for approximately 2 hour run time
    - 7. Programmable set point range of 65-180°F (18-82°C)
  - D. MEGATRON-PNV-125-LF shall be certified to ASSE standard 1017 and CSA B125.3-18 and so certified and identified
  - E. MEGATRON-PNV-125-LF shall be cULus listed and identified
  - F. Maximum Operating Pressure shall be no greater than 125 PSI
  - G. Disinfection Option with Proton Box 2.0
    - a. Optional 3T (3 additional temperature probes for hot and cold inlet in addition to return)
    - b. Optional REL (Connectivity for 5 relay states)
  - H. BMS connectivity available with Proton 2.5 box utilizing BACnet MS/TP connection in addition to including Disinfection
    - a. Optional 3T (3 additional temperature probes for hot and cold inlet in addition to return)
    - b. Optional REL (Connectivity for 5 relay states)
  - I. Wi-Fi Connectivity available with Proton 3.0 box in addition to BACnet MS/TP connection, Disinfection, relay, and 3 additional temperature probes for hot and cold inlet in addition to return
  - J. Megatron PNV-125-LF shall be strut mounted and include: outlet test connection with full port ball valve, Taco 0034ePLus circulator, boiler drain tee with analog temperature gauge, 1" recirculation loop tie-in with check valves and full port ball valves, 1" check valve on the cold and hot water Inlet with inlet analog temperature gauge and full port ball valve.